

Fulfilling your needs in the most demanding situations

Detection, Enumeration, Identification and Antimicrobial Susceptibility Testing

- Complete diagnosis in an all-in-1 product
- Compatible for use with liquid-Amies collection and transport media
- Optimized management of co-infections
- Antibiotics compliant with worldwide recommendations
- The right pathologic threshold tested for each specimen type
- Fast results for clinicians
- Visual, colorimetric reading



Complete diagnosis in 24h*

- Detection, Enumeration
- Identification
- Antimicrobial Susceptibility Testing

Cutting-edge performance

Sensitivity

- Detection and enumeration of *U. urealyticum* at 10³, 10⁴ and ≥10⁵ CCU/mL
- Detection and enumeration of M. hominis at 10⁴ CCU/mL
- The right pathologic threshold tested for each specimen type

Specificity

• Selective growth media for each species tested separately

Reliability

- 5 antibiotics (LVX/MXF/ERY/CLI/TET) tested in compliance with CLSI guidelines (Clinical & Laboratory Standards Institute)
- Doxycycline added to the panel

Fast results

• Integrated growth activator for M. hominis


Ease of use

- Ready-to-use reagents
- Each well pre-filled with growth medium
- · Visual colorimetric reading
- No instrument required

Flexibility and safety

- Optimized transport medium: Specimens preserved up to 72 hours before analysis
- Clip-on closing system for the tray

METHODOLOGY



^{*} Reading of all samples in 24h. Second reading needed for liquid samples having negative results within 24h

PRODUCT RANGE

	REF.	NAME	QTY	PURPOSE
	08000	MYCOFAST® RevolutioN 2	25 tests	Complete kit
Optional——	00081	MYCOFAST® RevolutioN 2 AMIES	25 tests	Complete kit
	00061	UMMt RevolutioN	50 vials	Additional transport medium
	00083	UMMt AMIES RevolutioN	50 vials	Additional transport medium (for Amies specimens)
	00082	COMPLEMENT MYCOFAST® RevolutioN 2	25 tests	Additional test trays
	00064	MYCOPLASMA Stabilizer	100 tests	Additive for sample storage > 56 hours
	00900	MYCOPLASMA Control	12 tests	U.u. and M.h. control strains

Product not available in the U.S Please contact your sales representative for product availability in your country





 Headquarters
 T: +33 1 41 45 07 10

 Australia
 T: +61 1800 815 098

 Belgium & Luxembourg
 T: +32 92 820 531

 Brazil
 T: +55 27 3025 1415

 France
 T: +33 4 83 36 10 82

 Italy
 T: +39 02 48 40 35 42

 New Zealand
 T: +64 0800 555 611

 Serbia
 T: +381 11 2467119

 Switzerland
 T: +41 26 663 86 60

 The Netherlands
 T: +31 313 430 581

 UK
 T: +44 1442 869320

 United States
 T: +1 609 216 7361

