

# MYCOPLASMA CONTROL

Ref. 00900



- 6 characterized strains *Ureaplasma urealyticum* and *Mycoplasma hominis* (Up to 12 tests)
- Long shelf life/ Easy-to-use
- Meets your quality management requirements (compliant with ISO 15189:2007 recommendations)
- Relevant results complying with M.I.C. values (Minimum Inhibitory Concentration) in liquid medium
- Identification of strains performed on U9 and Hayflick media
- Compatible use with MYCOFAST® EvolutionN 2, MYCOFAST® EvolutionN 3, and with main other kits in the market

## SIMPLE PROTOCOL

- 1 Reconstitute the U.u or M.h. strain with 1mL of sterile distilled water
- 2 Dilute the bacterial suspension into the transport medium or culture medium
- 3 Regenerate the growth medium if required
- 4 Inoculate the tray or microplate

PART NR: EMR 00900-490-2011/00EN



### WORLDWIDE OFFICES:

Belgium T: + 32 9 2820 531  
Brazil T: + 55 27 3025 1415  
France T: + 33 4 90 17 54 50  
Italy T: + 39 02 48 40 35 42

Middle East & Africa T: + 971 4 375 2744  
Serbia T: + 381 11 2467119  
Switzerland T: + 41 26 663 86 60  
The Netherlands T: + 31 485 511 675

UK T: + 44 1268 493377  
United States T: + 1 609 216 7360  
855-ELITECH (855-354-8324)

Headquarters T: + 33 1 41 45 07 13  
www.elitechgroup.com  
info@elitechgroup.com