

# Panoscreen<sup>®</sup>

Flexible



**Antibody Screening Cells** with flexible shipping schedules to fit your laboratory needs.

**IMMUCOR**

Transfuse | Transplant | Transform **a life**



# Panoscreen®

For the detection of unexpected red blood cell antibodies.

**Flexible shipping schedules based on your inventory needs**

- Manufactured weekly using fresh donors enables inventory of multiple in-date lots with unique cells
- Each cell vial contains a 2-4% suspension of single donor, group O red blood cells
- Industry exclusive primary lot of Panoscreen® I & II provides a homozygous expression of Fy<sup>a</sup> and Jk<sup>a</sup> on a 2-cell antibody screen
- All donor cells are screened for the Bg<sup>a</sup> antigen
- Hemantigen® – a convenient source of control cells for typing reagents

Panoscreen® I & II	Product Code:	0002380
	Size:	2 x 10 mL
	Specifications:	2-cell Antibody Screen R <sub>1</sub> R <sub>1</sub> and R <sub>2</sub> R <sub>2</sub> red cell donor phenotypes
Panoscreen® I & II, Primary Lot <i>Industry Exclusive</i>	Product Code:	0002390
	Size:	2 x 10 mL
	Specifications:	Primary lot of 2-cell Antibody Screen R <sub>1</sub> R <sub>1</sub> and R <sub>2</sub> R <sub>2</sub> red cell donor phenotypes Homozygous expression of Fy <sup>a</sup> , Jk <sup>a</sup>
Panoscreen® I, II and III	Product Code:	0002381
	Size:	3 x 10 mL
	Specifications:	3-cell Antibody Screen R <sub>1</sub> R <sub>1</sub> , R <sub>2</sub> R <sub>2</sub> and rr red cell donor phenotypes Homozygous expression of Fy, Jk, M, N, S, s
Hemantigen®	Product Code:	0002223
	Size:	1 x 10 mL
	Specifications:	Pooled Antibody Screen

To learn more or to set up a consultation with one of our Blood Bank Business Managers, call **855.IMMUCOR** (855.466.8267) or visit **www.immucor.com**.

**In Canada:**

t 1.800.565.0653

e DBLCustomerService@immucor.com

